












## Legacy Items: For use with BOA® dials

SKU	Product	Description	Picture	Available for warranty purposes only until
PK0060-000-10	RevoFit® Re-Lacing Kit	Contains 2m of legacy Lace compatible only with plastic BOA dials.		7/1/26
RP0100-320-05	RevoFit® Dial and Lace Replacement Kit	Replacement BOA dial and 2m of legacy Lace compatible only for legacy RevoFit systems.		7/1/26
LOK014	RevoFit® High Power Soft Lace - 10m Spool	Bulk replacement legacy High Power Lace (10m spool) compatible only with BOA dials and legacy RevoFit Kits. <b>NOT COMPATIBLE with Click® Reel</b>	N/A	7/1/26
PK3001	RevoLock® Dummy	Replacement Lamination Dummy compatible only with legacy RevoLock Kits using the BOA M3 dial.		1/1/22
RP3000-220-05	RevoLock® Lanyard Dial and Lace Replacement Kit	Replacement BOA dial and 1.2m of legacy Lace compatible only with legacy RevoLock Kits using the BOA M3 dial.		7/1/26

RP3001-000-05	RevoLock® Lanyard Re-Lacing Kit	Contains 1.2m of legacy Lace compatible with RevoLock sockets built with BOA dials.		7/1/26
RP3002-000-05	RevoLock® Lanyard - Lower-Ex Threaded Insert	Replacement RevoLock Lower-Extremity Threaded Insert which attaches the lace to the liner via a threaded connection. <b>NOT COMPATIBLE with RevoLock 4-Hole or RevoLock Align Kits.</b>		7/1/26
RP3003-000-05	RevoLock® Lanyard - Upper-Ex Threaded Insert	Replacement RevoLock Upper-Extremity Threaded Insert which attaches the lace to the liner via a threaded connection. <b>NOT COMPATIBLE with RevoLock 4-Hole or RevoLock Align Kits.</b>		7/1/26
LOK020	RevoLock® Soft Lace - 10m Spool	Bulk replacement legacy Spectra Lace (10m spool) compatible only with BOA dials and legacy RevoLock Kits. <b>NOT COMPATIBLE with Click® Reel</b>	N/A	7/1/26
RP0060-000-05	RevoFlex® Re-Lacing Kit	Replacement legacy Lace (2m) compatible only with BOA dials mounted on plastic hinged AFO's.		7/1/26

RP1100-320-05	RevoFlex® Dial and Lace Replacement Kit	Replacement BOA dial and 2m of legacy Lace compatible only with BOA dials mounted on plastic hinged AFO's.		7/1/26
LOK013	RevoFlex® Lace - 10m Spool	Bulk replacement legacy Lace (10m spool) compatible only with the BOA dials and legacy RevoFlex Kits.	N/A	7/1/26
RP0200-130-05	QuickFit® XS/SM/MD Dial and Lace Replacement Kit	Replacement BOA dial and 1.5m of compatible legacy Lace for XS/SM/MD size QuickFit Lacers using BOA dials.		7/1/26
RP0200-230-05	QuickFit® LG Dial and Lace Replacement Kit	Replacement BOA dial and 2m of compatible legacy Lace for LG size QuickFit Lacers using BOA dials.		7/1/26
LOK019	Low-Mid Power Steel Lace - 10m Spool	Bulk replacement legacy Low-Mid Power Steel Lace (10m spool) compatible only with BOA dials and legacy QuickFit Kits. <b>NOT COMPATIBLE with Click® Reel</b>	N/A	7/1/26