











## Legacy Products (Made with the BOA® dial)



# RevoFit®

SKU	Product	Description	Picture
PK2000-320-05	RevoFit Lamination Kit	With RevoFit, you can build in user adjustability almost anywhere on a device. It's adjustable through clothing and allows for innovative device design.	
PK2150-320-05	RevoFit Diagnostic Kit	RevoFit Diagnostic Kit lets you adjust even earlier in socket development. Great for patient buy-in, determining dial location and being confident in panel placement.	
PK2010-320-05	RevoFit FABRICATORS 10-pack	Large practices, regular users and central fabricators save money and always have a kit ready to use by ordering the Fabricators 10-Pack.	<ul style="list-style-type: none"> <li>• 10 Dials</li> <li>• 10 Lamination Dummies</li> <li>• 10 - 6ft sections of tubing/sheathing</li> <li>• 10 RevoFit Lacing Kits</li> <li>• 20 Meters RevoFit Lace</li> </ul>

RP2000-000-05	RevoFit Dummy w/ Silicone Overmold	Replacement Lamination Dummy	
PK0080-000-10	Bulk Tubing with Sheathing (6')	6' tubing for RevoFit	
PK0060-000-10	RevoFit® Re-Lacing Kit	Contains 2m of legacy Lace compatible only with plastic BOA dials.	
RP0100-320-05	RevoFit® Dial and Lace Replacement Kit	Replacement BOA dial and 2m of legacy Lace compatible only for legacy RevoFit systems.	
LOK014	RevoFit® High Power Soft Lace - 10m Spool	Bulk replacement legacy High Power Lace (10m spool) compatible only with BOA dials and legacy RevoFit Kits. <b>NOT COMPATIBLE with Click® Reel</b>	N/A



# RevoLock®

SKU	Product	Description	Picture
PK3000-220-05	RevoLock Lanyard	With RevoLock, you can build in user adjustable suspension with a socket. Offers automatic alignment, makes it easy to don while seated and reduces pistoning.	 A collection of RevoLock Lanyard components including two black BOA dials, a white lanyard, a metal ring, a small metal pin, and a metal rod.
PK3001	RevoLock Dummy	RevoLock Dummy comes with baseplate, lamination dummy and dummy screws.	 A black circular baseplate with two screws.
RP3000-220-05	RevoLock® Lanyard Dial and Lace Replacement Kit	Replacement BOA dial and 1.2m of legacy Lace compatible only with legacy RevoLock Kits using the BOA M3 dial.	 A collection of components including two black BOA dials, a white lanyard, a metal ring, a small metal pin, and a metal rod.
RP3001-000-05	RevoLock® Lanyard Re-Lacing Kit	Contains 1.2m of legacy Lace compatible with RevoLock sockets built with BOA dials.	 A collection of components including a white lanyard, a metal ring, and a small metal pin.



RP3002-000-05	RevoLock® Lanyard - Lower-Ex Threaded Insert (10mm)	Replacement RevoLock Lower-Extremity Threaded Insert which attaches the lace to the liner via a threaded connection. <b>NOT COMPATIBLE with RevoLock 4-Hole or RevoLock Align Kits.</b>	
RP3003-000-05	RevoLock® Lanyard - Upper-Ex Threaded Insert (6mm)	Replacement RevoLock Upper-Extremity Threaded Insert which attaches the lace to the liner via a threaded connection. <b>NOT COMPATIBLE with RevoLock 4-Hole or RevoLock Align Kits.</b>	
LOK020	RevoLock® Soft Lace - 10m Spool	Bulk replacement legacy Spectra Lace (10m spool) compatible only with BOA dials and legacy RevoLock Kits. <b>NOT COMPATIBLE with Click® Reel</b>	N/A









# RevoFlex

SKU	Product	Description	Picture
RP0060-000-05	RevoFlex® Re-Lacing Kit	Replacement legacy Lace (2m) compatible only with BOA dials mounted on plastic hinged AFO's. (Available for maintenance and warranty purposes only)	
RP1100-320-05	RevoFlex® Dial and Lace Replacement Kit	Replacement BOA dial and 2m of legacy Lace compatible only with BOA dials mounted on plastic hinged AFO's. (Available for maintenance and warranty purposes only)	
LOK013	RevoFlex® Lace - 10m Spool	Bulk replacement legacy Lace (10m spool) compatible only with the BOA dials and legacy RevoFlex Kits. (Available for maintenance and warranty purposes only)	N/A

# QuickFit

SKU	Product	Description	Picture
OK0180-130-05	QuickFit Lacer, XS	The extra small size QuickFit Lacer is designed for custom pediatric footwear size 35 EU and smaller. It comes with a low-profile dial and 4 guides. Rivets sold separately, but are required to install the QuickFit Lacer Kits. Lacers are sold individually. For a pair of shoes, two kits are required.	
OK0185-130-05	QuickFit Lacer, XS 10-Pack	The extra small size QuickFit Lacer is designed for custom pediatric footwear size 35 EU and smaller. It comes with 10 Extra Small QuickFit Lacer kits and 100 rivets. For a pair of shoes, two kits are required.	N/A
OK0200-130-05	QuickFit Lacer, SM	The small size QuickFit Lacer is designed for footwear size 36 - 41 EU. It comes with a low profile dial and 6 guides. Rivets sold separately, but are required to install the QuickFit Lacer Kits. Lacers are sold individually. For a pair of shoes, two kits are required.	

OK0205-130-05	QuickFit Lacer, SM 10-Pack	The small size QuickFit Lacer is designed for footwear size 36 - 41 EU. It comes with 10 Small QuickFit Lacer kits and 100 rivets. For a pair of shoes, two kits are required.	N/A
OK0220-130-05	QuickFit Lacer, MD	The medium size QuickFit Lacer is designed for footwear size 42 EU and larger. It comes with a low profile dial and 8 guides. Rivets sold separately, but are required to install the QuickFit Lacer Kits. Lacers are sold individually. For a pair of shoes, two kits are required.	
OK0225-130-05	QuickFit Lacer, MD 10-Pack	The Medium size QuickFit Lacer is designed for custom pediatric footwear size 42 EU and larger. It comes with 10 Medium QuickFit Lacer kits and 100 rivets. For a pair of shoes, two kits are required.	N/A
TL001	QuickFit Tool	The QuickFit Tool is customized for use with QuickFit Rivets and is the recommended installation tool. The QuickFit Tool properly secures rivets and provides either an 1/8" and 3/16" rivet hole.	
SP1040	QuickFit Rivets 100 Pack	QuickFit Rivets come in packs of 100. You can supply your own rivets, or purchase our easy to use rivet pack to install QuickFit Lacers to any device or shoe.	

RP0200-130-05	QuickFit® XS/SM/MD Dial and Lace Replacement Kit	Replacement BOA dial and 1.5m of compatible legacy Lace for XS/SM/MD size QuickFit Lacers using BOA dials.	
RP0200-230-05	QuickFit® LG Dial and Lace Replacement Kit	Replacement BOA dial and 2m of compatible legacy Lace for LG size QuickFit Lacers using BOA dials. (Available for maintenance and warranty purposes only)	
LOK019	Low-Mid Power Steel Lace - 10m Spool	Bulk replacement legacy Low-Mid Power Steel Lace (10m spool) compatible only with BOA dials and legacy QuickFit Kits. <b>NOT COMPATIBLE with Click® Reel</b>	N/A
OK1410-130-05	QuickFit Strap 1.5"	<ul style="list-style-type: none"> <li>• Prefabricated, durable, and install in minutes</li> <li>• Velcro® closure</li> <li>• 100mm (4") of adjustable travel</li> <li>• Usable strap length: Max: 0.75m (30") / Min: 0.2m (9")</li> <li>• Cut to size</li> <li>• Width (38mm)</li> </ul>	
OK1415-000-05	QuickFit Buckle 1.5"	<ul style="list-style-type: none"> <li>• Easily adds onto QuickFit Strap in seconds</li> <li>• Adds 100mm (4") to strap usable length</li> <li>• 5mm (.2") profile</li> <li>• Width (38mm)</li> </ul>	