

## With the **RevoFit®**, you can build in user adjustability *anywhere* on a socket.

### Micro/Macro Adjustments

Each "click" of the dial provides 1mm of movement for targeted compression. When used correctly, the system can adjust for up to 10 ply of sock changes.

### Adjustable Through Clothing

Operate the dial through clothing without removing or adding layers to adjust fit of the socket.



### Rapid Donning & Doffing

Releasing the tension in the system makes donning and doffing easy for bulbous ends or large volume fluctuations.

### Seated Comfort

Easily release the dial through clothing to make sitting more comfortable.

**Need user-friendly suspension solutions? Order RevoLock® Lanyard (PK3000-220-05).**

#### RevoFit® Diagnostic Kit: SKU #PK2150-320-05



#### RevoFit® Lamination Kit: SKU #PK2000-320-05



#### Lamination 10-Pack: SKU #PK2010-320-05





Save money and have the right kit in stock by ordering a 10-pack of RevoFit® Lamination Kits.

For additional tech specs, please go to: <https://clickmedical.co/legacy-products/>

# Replacement Parts

For additional tech specs, please go to:  
<https://clickmedical.co/legacy-products/>

A full line of replacement parts allows the system to remain functioning like new throughout the lifetime of the socket.

SKU	Item No.	Details
<b>RevoFit® Fabricators 10-pack</b> Large practices, regular users and central fabricators save money and always have a kit ready to use by ordering the Fabricators 10-Pack.	PK2010-320-05	Includes: <ul style="list-style-type: none"> <li>• 10 dials</li> <li>• 10 Lamination Dummies</li> <li>• 10 - 6ft sections of tubing/sheathing</li> <li>• 10 RevoFit Lacing Kits</li> <li>• 20 Meters RevoFit Lace</li> </ul>
<b>RevoFit® Dummy w/ Silicone Overmold</b> 	RP2000-000-05	Includes: <ul style="list-style-type: none"> <li>• Replacement Lamination Dummy</li> </ul>
<b>Bulk Tubing with Sheathing (6')</b> 	PK0080-000-10	Includes: <ul style="list-style-type: none"> <li>• 6' tubing for RevoFit</li> </ul>
<b>RevoFit® Re-Lacing Kit</b> 	PK0060-000-10	Includes: <ul style="list-style-type: none"> <li>• 7ft (2.1m) of 500-lb RevoFit® Lace</li> <li>• Lace Feeder Tool</li> <li>• T6 Tool</li> </ul> <p>We recommend replacing lace every 6 months or at any sign of abrasion to maintain proper system function.</p>
<b>RevoFit® Dial and Lace Replacement Kit</b> 	RP0100-320-05	Includes: <ul style="list-style-type: none"> <li>• RevoFit® Dial</li> <li>• 7ft (2.1m) of 500-lb RevoFit® Lace</li> <li>• Lace Feeder Tool</li> <li>• T6 Tool</li> </ul> <p>We recommend replacing lace at time of dial replacement.</p>
<b>RevoFit® High Power Soft Lace - 10m Spool</b>	LOK014	Includes: <ul style="list-style-type: none"> <li>• Bulk replacement legacy High Power 500-lb Lace (10m spool) compatible only with BOA dials and legacy RevoFit Kits.</li> </ul>