

With QuickFit[®] Straps and QuickFit[®] Buckles, you can easily add on micro-adjustable straps with a buckle that never wears out.

The QuickFit Strap

Micro-Adjustable

Empower patients to quickly and easily control fit and comfort throughout the day. Each "click" of the dial moves the strap 1mm for the perfect fit every time.

Improved Compliance

Patients LOVE having adjustability to create proper fit, which means they wear their devices more often.

One-Handed Use

Adjustment and closure can be operated with only one hand.



The QuickFit Buckle

Self-Locating Buckle

Powerful magnetic buckle clicks itself into place and holds strong all day.

One-Handed Use

Easily operate with only one hand.

Never Wears Out

The QuickFit Buckle design allows straps to last the life of the product by eliminating repetitive hooking and unhooking of Velcro[®].

| SKU | Size | Item No. | Details |
|---|-------------|---------------|--|
| QuickFit[™] Strap  | 1.5" / 38mm | OK1410-130-05 | <ul style="list-style-type: none"> • Prefabricated, durable, and install in minutes • Velcro[®] closure • 100mm (4") of adjustable travel • Usable strap length: <ul style="list-style-type: none"> • Max: 0.75m (30") / Min: 0.2m (9") • Cut to size • Made in USA |
| QuickFit[™] Buckle (sold separately)  | 1.5" / 38mm | OK1415-000-05 | <ul style="list-style-type: none"> • Easily adds onto QuickFit Strap in seconds • Adds 100mm (4") to strap usable length • Cut to size • 5mm (.2") profile • 70kg (150 lbs) holding capacity • Made in USA |