







Click Medical[®] Product Catalog

Click Medical products are patent protected. For full list see www.clickmedical.co/patents








SKU	Product	Description	Picture
RF-200-07-01	RevoFit® Lamination Kit with Click® Reel	<p>This single RevoFit® Lamination Kit with the Click® Reel allows you to fabricate user-controlled adjustability anywhere on a device.</p> <p>The Click® Reel innovation creates adjustable devices with a lower streamlined profile, micro/macro adjustment that incrementally tightens and loosens with increased durability.</p> <p>Kit includes: 1x Click Reel, 1x Lamination Collar, 2m HD Lace, 1x Metal Lace Feeder, 3x Plastic Lace Feeders, 1x Reel Tool, 1.8m Lamination Tubing</p> <p>User Weight Limit: 300lb (135kg)</p>	
RF-200-07-10	RevoFit® Lamination Fabricators 10-Pack with Click® Reel	<p>The RevoFit® Lamination Fabricators 10-Pack with the Click® Reel is recommended for central fabricators and practices that create large volumes of adjustable devices.</p> <p>Save 10% on purchase price of each single kit while assuring you will always have a kit on hand.</p>	
RP-002	HD Lace Kit	<p>2m of HD (High Durability) Lace and Click® Reel components for RevoFit® Lamination Kits with the Click® Reel.</p>	
RP-002-07	Click® Reel and HD Lace Kit	<p>Replacement Click® Reel and 2m of HD (High Durability) Lace for RevoFit® Kits only.</p> <p>Note: This is not for replacing a BOA® dial. Please see RP0100-320-05 for replacing a BOA® dial.</p>	

SKU	Product	Description	Picture
RP-200	Lamination Collar	<p>Replacement Lamination Collar for RevoFit® & RevoLock Lamination Kits with the Click® Reel only.</p>	
RP-201	Lamination Tube - 1.8m	<p>Replacement Lamination Tube (1.8m) for all Revo Kits with the Click Reel and BOA® dials.</p>	
BL-002	HD Lace - 25m Spool	<p>25m spool of HD (High Durability) Lace assures you will always have replacement lace on hand.</p> <p>Cost effective and sustainable with less waste.</p> <p>Compatible with both the RevoFit® Lamination Kit with the Click Reel and RevoFit® Lamination Kit with BOA® dial.</p>	
RP-300	Slider Part Kit	<p>Replacement Slider and Lock Wire for dynamic AFOs, large gap prostheses or RevoLock® applications.</p>	
RP-309	Slider Lamination Kit	<p>Lamination Kit for fabricating new dynamic AFOs, large gap prostheses or replacement for RevoLock® Kits.</p>	
AC-100-05	Comfort Grip Ring - Black 5-Pack	<p>Provide a soft-touch feel and 10% increased turning power to the existing Click® Reel. Great for patients with dexterity issues.</p>	
AC-101-05	Power Grip Ring - Black 5-Pack	<p>Provide a soft-touch feel and 20% increased turning power to the existing Click® Reel. Great for patients with dexterity and hand strength issues.</p>	
AC-000-00	Click® Care Kit	<p>The Click Care Kit enables you to service your clients in the office or on the road with everything readily accessible for maintaining, repairing or replacing Click components.</p> <p>Contents: 1x Surface Collar Kit, 2x HD Lace Kits, 2x HC Lace Kits, Flathead Screwdriver, Dental Pick, Spectra Scissors, Needle Nose Pliers & Reel Lubricant. Extra space is allocated for a Click® Reels and a lighter, which can be purchased separately.</p>	

RevoSurface®










SKU	Product	Description	Picture
RS-100-07-01	RevoSurface® Tool Kit with Click® Reel	<p>"This single RevoSurface Tool Kit with the Click Reel allows you to test adjustable designs and fit of check sockets, 3D printed devices, and on orthoses. It can also be used to solve fit issues on existing rigid devices.</p> <p>The Click Reel innovation creates adjustable devices with a lower streamlined profile, micro/macro adjustment that incrementally tightens and loosens with increased durability.</p> <p>Kit includes: 1x Click Reel, 2x Surface Collar, 2m HC Lace, 1.8m Lamination Tubing, 1x Metal Lace Feeder, 3x Plastic Lace Feeders, 1x Reel Tool, 6x Surface Guides, 3x Releasing Guides, 4x Chicago Screw assemblies</p> <p>User Weight Limit: 300lb (135kg)"</p>	
RP-001	HC Lace Kit	2m of HC (High Capacity) Lace and Click Reel components for RevoSurface Tool Kit, RevoLock Align, and RevoLock 4-Hole Kits with the Click Reel.	
RP-000-07	Click® Reel and HC Lace Kit	<p>Replacement Click Reel and 2m of HC (High Capacity) Lace for RevoLock Align, RevoLock 4-Hole and RevoSurface Tool Kits only.</p> <p>This is not for replacing a BOA® dial. Please see RP0100-320-05 for replacing a BOA® dial.</p>	
RP-101	Surface Guide 6-Pack	The 6-Pack of Surface Guides allows a practice to have on hand additional guides for easy replacement and configuration changes. For RevoSurface Tool Kits with Click Reel only.	





RevoSurface®


SKU	Product	Description	Picture
RP-102	Surface Collar	2x Click Reel mounting collar. For surface mount only. Not for use with lamination.	
RP-105	Surface Releasing Tab 3-Pack	3-Pack Surface Releasing Tabs for converting a device from a fixed lace path to a releasing lace path to fully open the device or to replace existing Surface Releasing Tabs.	
RP-201	Lamination Tube - 1.8m	Replacement Lamination Tube (1.8m) for all Revo Kits with the Click Reel and BOA® dials.	
BL-001	HC Lace - 25m Spool	<p>25m spool of HC (High Capacity) Lace assures you will always have replacement lace on hand.</p> <p>Cost effective, sustainable and less waste.</p> <p>Used in the RevoSurface, RevoLock Align, RevoLock 4-Hole, and RevoLock Lanyard V2 & V3 with BOA Dial.</p>	
AC-100-05	Comfort Grip Ring - Black 5-Pack	<p>Provide a soft-touch feel and 10% increased turning power to the existing Click Reel.</p> <p>Great for patients with dexterity issues.</p>	
AC-101-05	Power Grip Ring - Black 5-Pack	<p>Provide a soft-touch feel and 20% increased turning power to the existing Click Reel.</p> <p>Great for patients with dexterity and hand strength issues.</p>	
AC-000-00	Click® Care Kit	<p>The Click Care Kit enables you to service your clients in the office or on the road with everything readily accessible for maintaining, repairing or replacing Click components.</p> <p>Contents: 1x Surface Collar Kit, 2x HD Lace Kits, 2x HC Lace Kits, Flathead Screwdriver, Dental Pick, Spectra Scissors, Needle Nose Pliers & Reel Lubricant. Extra space is allocated for a Click® Reels and a lighter, which can be purchased separately.</p>	

SKU	Product	Description	Picture
RL-301-07-01	RevoLock® Align Kit with Click® Reel	<p>This single RevoLock Align Kit with the Click Reel is for transferable alignment fabrication for both upper and lower extremity sockets.</p> <p>The Click Reel innovation offers increased durability and power, a lower streamlined profile, and micro/macro adjustment that incrementally tightens and loosens. The Slider enhances safety for limbs changing in position or fluctuating in volume.</p> <p>Kit includes: 1x Click Reel, 1x Lamination Collar, 1x Slider Lamination Kit, 2x RevoLock LE Quick Connect Inserts, 1x RevoLock LE Quick Connect Snap, 2m HC Lace, 1x Metal Lace Feeder, 3x Plastic Lace Feeders, 1x Reel Tool, 0.5m Lamination Tubing, 1x Align Housing + Fabrication Post</p> <p>User Weight Limit: 300lb (135kg)</p>	
RL-300-07-01	RevoLock® 4-Hole Kit with Click® Reel	<p>This single RevoLock 4-Hole Kit with the Click Reel is for standard alignment fabrication for lower extremity sockets and comes with a standard 25mm x 65mm distal housing and integrated 4-Hole plate.</p> <p>The Click Reel innovation offers increased durability and power, a lower streamlined profile and micro/macro adjustment that incrementally tightens and loosens. The Slider enhances safety for limbs changing in position or fluctuating in volume. RevoLock® Tool Form required for fabrication of</p> <p>RevoLock 4-Hole kit. Kit includes: 1x Click Reel, 1x Lamination Collar, 1x Slider Lamination Kit, 2x RevoLock LE Quick Connect Inserts, 1x RevoLock LE Quick Connect Snap, 2m HC Lace, 1x Metal Lace Feeder, 3x Plastic Lace Feeders, 1x Reel Tool, 0.5m Lamination Tubing, 1x 4-Hole Puck, 1x Tool Form Plug</p> <p>User Weight Limit: 300lb (135kg)</p>	
RL-310-00-01	RevoLock® Tool Form	<p>The 25mm x 65mm RevoLock Fabrication Tool Form is a required, reusable fabrication tool for RevoLock® 4-Hole Kits.</p>	

SKU	Product	Description	Picture
RL-302-22-01	RevoLock® Lanyard V3 with BOA® Dial	<p>This kit offers the push to engage and pull to release interface of the BOA® Dial in a RevoLock Lanyard solution. It includes a shorter and smaller diameter distal housing which accepts the Quick Connect Snap for easy connection to the liner.</p> <p>Kit includes: 1x BOA® Lanyard Dial, 1x Lamination Collar, 1x 0.5m Lamination Tubing, 2x LE Quick Connect Inserts, 1x Quick Connect Snap, 1x Align Housing + Fabrication Post, 2x 1.22m HC Lace, 1x Metal Lace Feeder</p> <p>User Weight Limit: 220lb (100kg)</p>	
RP-310	RevoLock® V2 & V3 BOA® Dial and HC Lace Kit (Port C)	<p>Replacement 1x BOA® Lanyard Dial, 1.22m of HC Lace (High Capacity), 1x Metal Lace Feeder. For use with Port C Lace exit.</p> <p>Compatible only with RevoLock Lanyard V2 and V3 Kits</p>	
RP-310A	RevoLock® V2 & V3 BOA® Dial and HC Lace Kit (Port A)	<p>Replacement 1x BOA® Lanyard Dial, 1.22m of HC Lace (High Capacity), 1x Metal Lace Feeder. For use with Port A Lace exit.</p> <p>Compatible only with RevoLock Lanyard V2 and V3 Kits</p>	
RP-312	RevoLock® HC Lace Kit	<p>1.22m of HC Lace (High Capacity), 1x Metal Lace Feeder.</p> <p>Compatible only with RevoLock Lanyard V2 and V3 Kits</p>	
RP-001	HC Lace Kit	<p>2m of HC (High Capacity) Lace and Click Reel components for RevoLock Align, RevoLock® 4-Hole, and RevoSurface Tool Kits with the Click® Reel.</p>	
RP-000-07	Click® Reel and HC Lace Kit	<p>Replacement Click® Reel and 2m of HC (High Capacity) Lace for RevoLock Align, RevoLock® 4-Hole and RevoSurface Tool Kits only.</p> <p>This is not for replacing a BOA® dial. Please see RP0100-320-05 for replacing a BOA® dial.</p>	

SKU	Product	Description	Picture
RP-200	Lamination Collar	Replacement Lamination Collar for RevoLock Align, RevoLock 4-Hole, RevoLock Lanyard V3 & RevoFit Lamination Kits.	
RP-201	Lamination Tube - 1.8	Replacement Lamination Tube (1.8m) for all Revo Kits with the Click Reel and BOA® dials.	
RP-305	RevoLock® Tool Form Plug 5-Pack	The RevoLock Tool Form Plug 5-Pack is to ensure fabricators have extra plugs on hand.	
RP-300	Slider Part Kit	Replacement Slider and Lock Wire for dynamic AFOs, large gap prostheses or RevoLock applications.	
RP-309	Slider Lamination Kit	Lamination Kit for fabricating new dynamic AFOs, large gap prostheses or replacement for RevoLock Kits.	
BL-001	HC Lace - 25m Spool	25m spool of HC (High Capacity) Lace assures you will always have replacement lace on hand Cost effective, sustainable and less waste. Used in the RevoSurface and RevoLock Kits with the Click Reel only.	
AC-100-05	Comfort Grip Ring - Black 5-Pack	Provide a soft-touch feel and 10% increased turning power to the existing Click Reel. Great for patients with dexterity issues.	
AC-101-05	Power Grip Ring - Black 5-Pack	Provide a soft-touch feel and 20% increased turning power to the existing Click Reel. Great for patients with dexterity and hand strength issues.	
RP-306	RevoLock® UE Fixed Insert	RevoLock Fixed Threaded Insert (M6) for upper extremity and pediatric liners. Compatible with RevoLock 4-Hole, RevoLock Align and RevoLock Lanyard V3 Kits. Best for patients who cannot operate the RevoLock Quick Connect Snap or for suction systems needing an airtight seal. This part is NOT releasable and does not connect to the RevoLock LE Quick Connect Snap.	

SKU	Product	Description	Picture
RP-307	RevoLock® LE Fixed Insert	RevoLock Fixed Threaded Insert (M10) for lower extremity liners only. Compatible with RevoLock 4-Hole, RevoLock Align and RevoLock Lanyard V3 Kits. Best for patients who cannot operate the RevoLock Quick Connect Snap or suction systems needing an airtight seal. This part is NOT releasable and does not connect to the RevoLock LE Quick Connect Snap.	
RP-303	RevoLock® LE Quick Connect Snap	Replacement RevoLock LE Quick Connect Snap connects the HC (High Capacity) Lace to the liner via the RevoLock LE Quick Connect Insert (RP-302).	
RP-302	RevoLock® LE Quick Connect Insert	Replacement RevoLock LE Quick Connect Insert for connection to lower extremity liners. Compatible with RevoLock LE Quick Connect Snap (RP-303) found in the RevoLock Align, RevoLock 4-Hole and RevoLock Lanyard V3 Kits. Often purchased for outfitting additional patient liners.	
AC-000-00	Click® Care Kit	The Click Care Kit enables you to service your clients in the office or on the road with everything readily accessible for maintaining, repairing or replacing Click components. Contents: 1x Surface Collar Kit, 2x HD Lace Kits, 2x HC Lace Kits, Flathead Screwdriver, Dental Pick, Spectra Scissors, Needle Nose Pliers & Reel Lubricant. Extra space is allocated for a Click® Reels and a lighter, which can be purchased separately.	








SKU	Product	Description	Picture
QFL-130-B-2	QuickFit® Lacer - Black 2-Pack	<p>The QuickFit Lacer is a dial-based lacing system that tightens a high-strength lace to provide adjustable closure for upper extremity prostheses, orthoses, and footwear. The low-profile dial enables smooth, one-handed micro-adjustments while distributing tension evenly across the closure area.</p> <p>Redesigned for easy integration, the system comes partially pre-assembled and allows practitioners to customize lace routing by placing guides where needed—no system disassembly required. QuickFit Lacers are sold in pairs, making them ideal for footwear and bilateral applications.</p> <p>Kit Includes: 2 × Low Power Dials with one end of 1 m black lace attached, 2 × Dial Housings, 20 × Black QuickLace Guides, 32 × Black Rivets</p>	
QFL-130-W-2	QuickFit® Lacer - White 2-Pack	<p>The QuickFit Lacer is a dial-based lacing system that tightens a high-strength lace to provide adjustable closure for upper extremity prostheses, orthoses, and footwear. The low-profile dial enables smooth, one-handed micro-adjustments while distributing tension evenly across the closure area.</p> <p>Redesigned for easy integration, the system comes partially pre-assembled and allows practitioners to customize lace routing by placing guides where needed—no system disassembly required. QuickFit Lacers are sold in pairs, making them ideal for footwear and bilateral applications.</p> <p>Kit Includes: 2 × Low Power Dials with one end of 1 m white lace attached, 2 × Dial Housings, 20 × White QuickLace Guides, 32 × White Rivets</p>	
QFL-230-B-2	QuickFit® Lacer Large - Black 2-Pack	<p>The QuickFit Lacer Large is a dial-based lacing system that tightens a high-strength lace to deliver powerful, adjustable closure for above-the-ankle footwear and orthoses. Designed specifically for applications such as leather AFOs, rigid plastic AFOs, and CROW boots, it provides dependable closure with a simple turn of a dial.</p> <p>Redesigned for easy integration, the system comes partially pre-assembled and allows practitioners to customize lace routing by placing guides where needed. QuickFit Lacers are sold in pairs, making them ideal for bilateral applications.</p> <p>Kit Includes: 2 × Mid Power Dial with a single end of 1.3 m of black lace attached, 2 × Dial Housings, 20 × Black QuickLace Guide, 32 × Black Rivets</p>	

SKU	Product	Description	Picture
RP-103-B	QuickLace Guide - Black 10-Pack	<p>QuickLace Guides are surface-mounted guides designed to route and manage the lace used with the QuickFit Lacer system. They help distribute tension across the closure area while allowing smooth, consistent lace movement.</p> <p>The guides attach easily to prostheses, orthoses, or footwear using rivets and can be positioned to customize the closure configuration. Their dual-function design allows the lace to be either permanently retained for secure routing or releasable so the device can fully open for easier donning and doffing.</p> <p>QuickLace Guides also provide a convenient option for replacing existing guides or adding additional lace crossings to adjust a device's closure performance.</p>	
RP-103-W	QuickLace Guide - White 10-Pack	<p>QuickLace Guides are surface-mounted guides designed to route and manage the lace used with the QuickFit Lacer system. They help distribute tension across the closure area while allowing smooth, consistent lace movement.</p> <p>The guides attach easily to prostheses, orthoses, or footwear using rivets and can be positioned to customize the closure configuration. Their dual-function design allows the lace to be either permanently retained for secure routing or releasable so the device can fully open for easier donning and doffing.</p> <p>QuickLace Guides also provide a convenient option for replacing existing guides or adding additional lace crossings to adjust a device's closure performance.</p>	
TL-001	QuickFit® Tool	<p>The QuickFit Tool is customized for use with QuickFit Rivets and is the recommended installation tool. The QuickFit Tool properly secures rivets and provides either an 1/8" and 3/16" rivet hole.</p>	
SP1040	QuickFit® Rivets - Black 100-Pack	<p>Color-matched rivets designed for securely attaching QuickLace Guides and dials when installing the QuickFit Lacer system. Ideal for new builds, replacements, or configuration changes across prosthetic, orthotic, and footwear applications.</p>	




QuickFit®

Made with the BOA® Dial





SKU	Product	Description	Picture
SP1040-S	QuickFit® Rivets - Silver 100-Pack	Color-matched rivets designed for securely attaching QuickLace Guides and dials when installing the QuickFit Lacer system. Ideal for new builds, replacements, or configuration changes across prosthetic, orthotic, and footwear applications.	
RP0200-130-05	QuickFit® Dial and Lace Replacement Kit - Black	Replacement dial and lace for the QuickFit Lacer system in black. The low-profile dial tightens the high-strength lace to provide smooth, precise adjustability for prostheses, orthoses, and footwear. Kit Includes: 1 × Low Power Dials with one end of 1 m black lace attached	
RP0200-130-W	QuickFit® Dial and Lace Replacement Kit - White	Replacement dial and lace for the QuickFit Lacer system in white. The low-profile dial tightens the high-strength lace to provide smooth, precise adjustability for prostheses, orthoses, and footwear. Kit Includes: 1 × Low Power Dials with one end of 1 m white lace attached	
RP0200-230-05	QuickFit® LG Dial and Lace Replacement Kit - Black	Replacement dial and lace for the QuickFit Lacer LG system in black. The mid-power dial tightens a high-strength lace to deliver powerful, adjustable closure for above-the-ankle footwear and orthoses. Kit Includes: 1 × Mid Power Dial with a single end of 1.3 m of black lace attached	
LOK019	Steel Lace - Black 10m Spool	The Steel Lace – Black 10 m Spool provides durable, color-matched lace for QuickFit Lacer systems, allowing practitioners to cut custom lengths while maintaining a clean, consistent system appearance.	

QuickFit®



Made with the BOA® Dial

SKU	Product	Description	Picture
LOK019-W	Steel Lace - White 10m Spool	The Steel Lace – White 10 m Spool provides durable, color-matched lace for QuickFit Lacer systems, allowing practitioners to cut custom lengths while maintaining a clean, consistent system appearance.	
OK1410-130-05	QuickFit® Strap 1.5"	QuickFit Straps are designed to quickly install on a wide variety of cuffs, orthoses, and prosthetic devices to create an improved, easy to use, adjustable strap. <ul style="list-style-type: none"> • Prefabricated, durable, and installs in minutes • Don/doff via simple Velcro® closure • 76mm (3") of adjustable travel • Usable strap length: Max: 0.75m (30") / Min: 0.2m (9") • Cut to size • Width (38mm) 	
OK1415-000-05	QuickFit® Buckle 1.5"	The QuickFit Buckle can be added onto the QuickFit Strap in a few seconds. The magnetic buckle allows the user to easily connect the strap - even with limited dexterity or in hard to reach areas. Its plastic on plastic interface eliminates wear issues of Velcro® closure. <ul style="list-style-type: none"> • Easily adds onto QuickFit Strap in seconds • Adds 76mm (3") to strap usable length • 5mm (.2") profile • Width (38mm) 	

Replacement Parts

SKU	Product	Description	Picture
PK0060-000-10	RevoFit® Re-Lacing Kit for BOA® dial	2m of HD (High Durability) Lace for RevoFit Lamination Kits. Compatible only with RevoFit Kits with BOA® dial.	
RP0100-320-05	RevoFit® Dial and Lace Replacement Kit with BOA® dial	Replacement BOA® dial and 2m of HD (High Durability) Lace Compatible only with RevoFit Kits with BOA® dial.	
RP3000-220-05	RevoLock® Lanyard Dial and Lace Replacement Kit with BOA® dial	Replacement BOA® dial and 1.2m of lace Compatible only with RevoLock Kits with BOA® dial. Not compatible with RevoLock V2 and V3 Kits.	
RP3001-000-05	RevoLock® Lanyard Re-Lacing Kit for BOA® dial	Contains 1.2m of lace Compatible only with RevoLock® Kits with BOA® dial. Not compatible with RevoLock® V2 and V3 Kits.	

Replacement Parts

SKU	Product	Description	Picture
RP3002-000-05	RevoLock® Lanyard - Lower-Ex Threaded Insert (10mm)	Replacement RevoLock Lower-Extremity Threaded Insert which attaches the lace to the liner via a threaded connection. Compatible only with RevoLock Lanyard and RevoLock Lanyard V2 Kits with BOA® dial.	
RP3003-000-05	RevoLock® Lanyard - Upper-Ex Threaded Insert (6mm)	Replacement RevoLock Upper-Extremity Threaded Insert which attaches the lace to the liner via a threaded connection. Compatible only with RevoLock Lanyard and RevoLock Lanyard V2 Kits with BOA® dial.	

Find Click Medical® Adjustable Solutions Worldwide



DISTRIBUTOR	COUNTRIES	ADDRESS	PHONE	WEB ADDRESS
NORTH AMERICA				
SPS	United States (includes Puerto Rico)	6530 Corporate Court Alpharetta, GA 30005 United States	Tel: +1 (800) 767-7776 ext.3	spsco.com
Cascade	United States (includes Puerto Rico) / Canada	2638 Aztec Drive Chico, CA 95928 United States	Tel: +1 (800) 888-0865	cascade-usa.com
OrtoPed	Canada	373 Rue McCaffrey Montréal, QC H4T 1Z7 Canada	Tel: +1 (514) 342-6161	ortoped.ca
PEL Supply	United States	4666 Manufacturing Avenue Cleveland, OH 44135-2672 United States	Tel: 1-800-321-1264	pelsupply.com
Ortiz International	Mexico	Federalismo Sur 836 Guadalajara, Jalisco, 44190 Mexico	Tel: +33 3811-2151 6 3811-3827	ortizinternacional.com.mx
EUROPE				
Alfimed AS	Norway, Sweden, Denmark, Finland	Olav Brunborgs vei 6 1396 Billingstad Norway	Tel: +47 31 29 68 86	alfimed.no
EMO	Spain, Portugal	Camino Pascualeta, 40. 46200 Paiporta, Spain	Tel: +34 963 132 051	emo.es/es
Frey Orthopädie-Bedarf AG	Switzerland	Panoramaweg 35 5504 Othmarsingen Switzerland	Tel: +41 62 887 45 00	freyortho.ch/en
Janton Solutions	France	Chemin des Aubrières 49400 Saumur France	Tel: +(33) 2 41 50 50 91	janton-g2m.com
ING Corporation	Czech Republic, Romania, Slovakia, Croatia, Latvia	Mánesova 1259, 739 11 Frýdlant nad Ostravicí Czech Republic	Tel: +(420) 558 663 217	ingcorporation.cz
Mediroyal	Norway, Sweden, Denmark, Finland, Iceland, Estonia	Box 7052 Staffans väg 6 B S-192 07 Sollentuna Stockholm, Sweden	Tel: +46 8 506 766 00	mediroyal.se/en

Metiz IMPEX	Russia	Kolpakova 24/A, Mytischy, Moscow region, 141008 Kolpakova 24/A Mytischy, Moscow region 141008 Russia	Tel: +7 (495) 215-53-82	metiz-ltd.com
BeMedico	Belgium, Netherlands, Luxembourg	Mutsaardstraat 47 9550 Herzele, Belgium	Tel: +32 (0) 54 50 40 05	bemedico.be
Ottobock	Germany, (worldwide)	Max-Näder-Straße 15 (or Max-Naeder-Strasse 15) D-37115 Duderstadt, Germany	Tel: +49 5527 848-0	ottobock.com
Protheseus GmbH	Germany, Austria	Liebigstraße 13 91126 Schwabach, Germany	Tel: +49 (0) 9122 88 707-0	protheseus.de
Invapol	Ukraine	Akademika Kyms'koho 27 Kyiv, UA, 3680	Tel: +38 0502321540	invapol.com
Simex	Ukraine	Building 3, Letter "U", str. Glubochitska 17 Kiev, Ukraine, 04050	Tel: +(38) 044 4861406	parasind.com
Steeper	UK, Ireland	Unit 3, Stourton Link Intermezzo Drive Leeds LS10 1BL UK	Tel: (+44) 0 113 270 4841	steepergroup.com
MIDDLE EAST				
Gilron	Israel	2 Haron Katzir Ave Tel-Hashomer, Ramat-Gan, Israel	Tel: +(972) 3-5409333	gilron.com
OrthoMENA	United Arab Emirates	Dubai Real Estate Centre 153 Al Mina road 1st Floor, 103 Al Hudaiba, Dubai	Tel: 0097148357207	orthomena.com
AUSTRALIA				
OttoBock Australia Pty. Ltd.	Australia, New Zealand	Suite 1.01, Century Corporate Centre 62 Norwest Boulevard, Baulkham Hills, NSW, 2153, Australia	Tel: +61 2 8818-2800	ottobock.com.au
ASIA				
Imasen Engineering Corporation	Japan	3-1-8 Techno Plaza Kakamigahara Gifu JAPAN 509-0109	Tel: +(81) 90-2348-9484	imasengiken.co.jp
Teh Lin Prosthetic & Orthopaedic Inc.	China, Taiwan	No.7, Wuquan 7th Rd., Wugu Dist., New Taipei City 248 New Taipei City, Taiwan 248 Taiwan, China	Tel: +(886) 2-22991901	tehlin.com
SOUTH AMERICA				
BM Medicine SRL	Argentina	Bernardo De Irigoyen 1312 Buenos Aires, Capital Federal 1138 Argentina	Tel: +5491136252594	bmmedicine.com.ar
Ethnos	Brazil	Rua 15 de Novembro, 469, Vila São Luiz, Duque de Caxias Rio de Janeiro, Brazil	Tel: +55 21 2775 4032	ethnos.com.br
Orthopedix Argentina	Argentina	Av Cordoba, 5328,1ºA, Buenos Aires 1414, Argentina	Tel: +541147758208	orthopedixargentina.com.ar

Innovation, Outcomes, Reimbursement

Reimbursement varies by country.

There are several International codes established. Please ask your distributor for information specific to your country.

RevoFit® and RevoSurface®

Custom Lower Extremity

Prosthesis - L5783: Addition to lower extremity, user adjustable, mechanical, limb volume management system. Applicable to all LE amputation levels across all activity levels. Appropriate for new, replacement, and preparatory sockets.



Custom Upper Extremity Prosthesis - New Code Effective

04/01/2025, L7406: Addition to upper extremity, user adjustable, mechanical, residual limb volume management system to describe RevoFit®. Prior to 04/01/2025, continue to utilize L7499 (Reference L5783 as Same/Similar and pending activation of L7406), and, supply Click Medical's RevoFit MSRP for upper extremity prosthetic devices and our suggested descriptive language.

Custom Orthosis - L2999 (Reference L5783 as Same/Similar):

Additionally, utilize Click Medical's RevoFit MSRP for orthotic devices and our suggested descriptive language. Appropriate for new, replacement, and modified Dynamic AFO and KAFO custom laminated devices.

Repair of an Existing Prosthesis - L7510: For cost of materials and L7520 for repair labor components. Appropriate for RevoFit replacement parts and labor to install.

Repair of an Existing Orthosis - L4210: For cost of materials and L4205 for repair labor components. Appropriate for RevoFit replacement parts and labor to install.

RevoLock® 4-Hole & Align

Lower Extremity Prosthesis - L5671: Addition to lower extremity, below knee/above knee suspension locking mechanism (shuttle, lanyard, or equal), excludes socket insert.

Upper Extremity Prosthesis - L6698: Addition to upper extremity prosthesis, below elbow/above elbow, lock mechanism, excludes socket insert.

Repair of an Existing Prosthesis - L7510: For cost of materials and L7520 for repair labor components. Appropriate for RevoLock replacement parts and labor to install.

Repair of an Existing Orthosis - L4210: For cost of materials and L4205 for repair labor components. Appropriate for RevoLock replacement parts and labor to install.

QuickFit®

Lower Extremity Custom Orthosis - L2320: Addition of lacer to lower extremity, non-molded lacer, for custom fabricated orthosis only.



Lower Extremity Custom Orthosis - L2330: Addition of lacer to lower extremity, lacer molded to patient model, for custom fabricated orthosis.

Orthotic Shoes - L2999 (Reference L3580 as Same/Similar): Addition of lacer to orthopedic shoe, convert instep to velcro closure.

Addition to/Repair of Orthosis - L4002: Replacement strap for any orthosis that includes all components and can be any length or type.

CMS Approves **L5783** & **L7406** for Click Medical

Expanding RevoFit® Coverage to ALL Prosthetics!

HCPCS CODE	COVERAGE	EFFECTIVE DATE	REIMBURSEMENT
L5783	 Lower Extremity	April, 1 2024	\$3,088.30
L7406	 Upper Extremity	April, 1 2025	\$3,088.30

Questions?

We are here for you!

Technical Support

Monday - Thursday, 9am - 5pm MST
Friday, 9 am - 3 pm MST
+1 (970) 670-7012

Online Support Form:

clickmedical.co/contact-us

Click Medical

1205 Hilltop Parkway, W101
PO Box 775245
Steamboat Springs, CO
80487



**Better for your patient.
Better for your practice.
Better for your **bottom line.****



clickmedical.co

