

Steps:

1. Choose the adjustable socket design
2. Indicate the location of the adjustable system
3. Select the reel location

Clinical Solutions for Hinge Design: Expansion for doff and don, targeted expansion

**Hinge Socket Design**

Check box for hinge design. (Hinge size and shape must also be indicated on check socket or socket file.) Hinge element: 3/16"x1"x1/2" flexible plastic unless otherwise indicated.




Hinged Open Posterior




Supracondylar

Alternative Hinge Design Specifications Notes: \_\_\_\_\_

**Reel Location**

Check box for reel location. (Reel placement must also be indicated on the check socket or socket file.)

Frame mounted:

<input type="checkbox"/>	or	<input type="checkbox"/>	or	<input type="checkbox"/>
Proximal		Mid-shaft		Distal
_____ & _____				
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Posterior	Anterior	Medial	Lateral	

**OR**

Hinged Panel Mounted

Alternative Reel Placement Notes: \_\_\_\_\_